

EXP7: Partial molar volume

1. Although to a great extent the weights themselves enter into the calculations, it is necessary also to have the density d of every solution to within an accuracy of at least one part per thousand:

$$d = W_{\text{sol}}/V = W - W_e/V_p$$

The volume of the pycnometer V_p is obtained by use of the density of pure water at 30.0°C , d_0 (with the value $0.9957 \text{ g} \cdot \text{cm}^{-3}$), and $W_0 - W_e$.

(sol) 先計算 M

(1). 由 H_2O 計算 $V_p = W - W_e/d = ?$ 其中: $W = W_{\text{sol}} + W_e$ W_e : 比重瓶重 d : 查表

(2). 分別計算 4 瓶 $d = W - W_e/V_p$

dilution	1/8	1/4	1/2	1
$M_{\text{NaCl}}(M)$				
$d(\text{g}/\text{cm}^3)$				

2. The molalities m (concentration in mole per Kg of solvent) which are needed for the calculations can be obtained from the molarities M (concentration in mole per liter of solution) obtained from the volumetric procedures by using the equation

$$m = 1/(M/d \cdot M_2/1000) \cdot M/d = 1/(d/M - M_2/1000)$$

Where M_2 is the solute molecular weight ($58.45 \text{ g}/\text{cm}$) and d is the experimental density in $\text{g} \cdot \text{cm}^{-3}$ units.

(sol) $M_2 = 58.45 \text{ g/mol}$

由 $m = 1/(d/M - M_2/1000)$ 計算出 4 個 $m(\text{mol}/\text{Kg})$

3. Calculate Φ for each solution using Eq. (7) and plot Φ vs. $m^{1/2}$. Determine the slope $d\Phi/dm^{1/2}$ and the intercept Φ^0 at m equal zero from the best straight line through these data points.

(sol)

(1) 計算出 4 個 $\Phi(\text{ml}/\text{mol})$

$$\Phi = 1/d (M_2 - (1000/m)(W - W_0)/(W_0 - W_e)) \quad \dots\dots\dots (7)$$

$$\text{或} = 1/d (M_2 - (1000/m)(d - d_0)/(d_0))$$

(2) Φ 對 $m^{1/2}$ 作圖

$\Phi(\text{ml}/\text{mol})$				
$m(\text{mol}/\text{kg})$				
$m^{1/2}(\text{mol}/\text{kg})^{1/2}$				

4. Calculate ψ_2 and ψ_1 for $m = 0, 0.5, 1.0, 1.5, 2.0$, and 2.5 . Plot them against m and draw a smooth curve for each of the two quantities

(sol)

(1)

$$\psi_1 = V_1^0 - m/55.51(m^{1/2}/2 \cdot d\Phi/dm^{1/2}) \dots\dots\dots(12)$$

$$\psi_2 = \Phi^0 + 3/2m^{1/2} \cdot d\Phi/dm^{1/2} \dots\dots\dots(11)$$

再由上式 Φ 對 $m^{1/2}$ 作圖得 $\Phi = a m^{1/2} + b$

可得 $m = 0$ 時 $\Phi_0 =$

$$d\Phi/dm^{1/2} =$$

(2) 由 Φ_0 、 $d\Phi/dm^{1/2}$ 、及 $m = 0, 0.5, 1.0, 1.5, 2.0$, and 2.5 . 分別代入(12)、(11)

求得 6 組 $\psi_1 \psi_2$

5. In your report, present the curves (Φ vs. m , $V^{1/2}$ and ψ_1 vs. m) mentioned above. Present also in tabular form the quantities d , M , m , and Φ for each solution studied. Give the values obtained for the pycnometer volume V_p and for Φ^0 and $d\Phi/dm^{1/2}$

dilution	H ₂ O	1/8	1/4	1/2	1
$M_{\text{NaCl}}(M)$					
$d(\text{g/ml})$					
$m(\text{mol/kg})$					
$\Phi(\text{ml/mol})$					

ψ_1 VS. m 作圖一張及 ψ_2 VS. m 作圖一張

$m(\text{mol/kg})$	0	0.5	1	1.5	2	2.5
$\psi_1(\text{cm}^3/\text{mol})$						

$m(\text{mol/kg})$	0	0.5	1	1.5	2	2.5
$\psi_2(\text{cm}^3/\text{mol})$						

Exp.7

Wt. of NaCl	34.998 g	room temperature	23.9°C
Molar conc. of NaCl	2.9944mole/L	room pressure	747.1 mmHg
W _e	27.197g		
T _{exp't}	23.8°C		
d _{water,exp't}	0.9974g/ml		

Run	H ₂ O	1/8	1/4	1/2	1
M _{NaCl}	0	0.3743	0.7486	1.4972	2.9944
W _{sol'n + We}	77.123	77.915	78.674	80.137	82.888

Calculation (以 3 M 為例)

1. NaCl_(aq) 取 34.998 g

$$M_{NaCl} = \frac{n_{NaCl}}{V(L)} = \frac{\frac{34.998(g)}{58.44\left(\frac{g}{mol}\right)}}{0.200(L)} = 2.994(M)$$

2. 比重瓶體積

$$V_p = \frac{W_{water}}{d_{water}} = \frac{(W_{sol'n} + W_e) - W_e}{d_{water}} = \frac{77.145 - 27.197}{0.9974} \\ = 50.078(ml)$$

3. 密度(density) (以 3M 為例)

$$d_1 = \frac{w_1}{V_p} = \frac{(82.888 - 27.197)g}{50.078 ml} = 1.1121 \left(\frac{g}{ml}\right)$$

4. Molalities (m) (以 3M 為例)

$$m_1 = \frac{1}{\frac{d}{M} - \frac{M_{NaCl}}{1000}} = \frac{1}{\left(\frac{1.1121\left(\frac{g}{ml}\right)}{2.9944\left(\frac{mol}{L}\right)}\right) - \left(\frac{58.44\left(\frac{g}{mol}\right)}{1000\left(\frac{g}{kg}\right)}\right)} = 3.195 \left(\frac{mol}{kg}\right)$$

5. ϕ (以 3M 為例)

$$\phi_1 = \frac{1}{d_1} \times \left(M_{NaCl} - \frac{1000}{m} \times \frac{d_1 - d_0}{d_0} \right)$$

$$\begin{aligned}
&= \frac{1}{1.1121\left(\frac{g}{ml}\right)} \times \left(58.44\left(\frac{g}{mol}\right) - \frac{1000}{3.195\left(\frac{mol}{kg}\right)}\right) \\
&\times \left(\frac{1.1121\left(\frac{g}{ml}\right) - 0.9974\left(\frac{g}{ml}\right)}{0.9974\left(\frac{g}{ml}\right)} \right) \\
&= 20.19_4\left(\frac{ml}{mol}\right)
\end{aligned}$$

6. ϕ vs. \sqrt{m} (y 為 ϕ x 為 \sqrt{m}) $y = 3.2448x + 14.237$


$$\frac{d\phi}{d\sqrt{m}} = slope = 3.2448 \quad m=0 \text{ 時 } \phi^0 = 14.237 \left(\frac{ml}{mol}\right)$$

7. V_1 & V_2 以 $m=0.5$ 為例

$$\begin{aligned}
\bar{V}_1 &= \bar{V}_1^0 - \frac{m}{55.51(mole)} \times \left(\frac{\sqrt{m}}{2} \times \frac{d\phi}{d\sqrt{m}} \right) \\
&= 18.063(ml) - \frac{0.5\left(\frac{mole}{kg}\right)}{55.51\left(\frac{mole}{kg}\right)} \\
&\times \left(\frac{\sqrt{0.5}\left(\frac{mol^{0.5}}{kg^{0.5}}\right)}{2} \times 3.2448\left(\frac{kg^{0.5}}{mol^{0.5}}\right) \right) \\
&= 18.053\left(\frac{ml}{mol}\right) \\
\bar{V}_2 &= \phi^0 + \left(\frac{3\sqrt{m}}{2} \times \frac{d\phi}{d\sqrt{m}} \right) \\
&= 14.237(ml) + \frac{3 \times \sqrt{0.5}\left(\frac{mole^{0.5}}{kg^{0.5}}\right)}{2} \times 3.2448\left(\frac{kg^{0.5}}{mol^{0.5}}\right) \\
&= 17.679\left(\frac{ml}{mol}\right)
\end{aligned}$$

$$(\bar{V}_1^0 = \frac{18.016}{d_{water}} = \frac{18.016}{0.9974} = 18.063 \left(\frac{cm^3}{mol}\right))$$

Run	H ₂ O	1/8	1/4	1/2	1
M _{NaCl}	0	0.3743	0.7486	1.4972	2.9944
W _{sol'n + We}	77.123	77.915	78.674	80.137	82.888
We			27.197		
W _{sol'n}	49.926	50.718	51.477	52.94	55.691
dwater (g/ml)			0.9974		
V _p (ml)			50.056		
d _{sol'n}	0.9974	1.0132	1.0284	1.0576	1.1126
m(mol/kg)	0	0.378	0.760	1.543	3.194
φ	0	16.21	17.09	18.27	20.03
√m	0.000	0.614	0.872	1.242	1.787



$d\psi/dm =$	3.2448		
$\psi^\circ =$	14.237	ml/mol	
m	V1	V2	\sqrt{m}
0.0	18.063	14.237	0.000
0.5	18.053	17.679	0.707
1.0	18.034	19.104	1.000
1.5	18.009	20.198	1.225
2.0	17.980	21.120	1.414
2.5	17.947	21.933	1.581

